



Buy Direct and Save | (800) 474-7665 | www.ExactaDental.com | info@exactadental.com



Why

EXACTA Dental was founded with a single purpose - to provide value to our customers.

Our products have been carefully designed or selected to provide you with outstanding performance and results. We maintain a direct relationship with you so that we can become familiar with your products, techniques, and preferences, so that ultimately we can serve you better.

For more info visit www.exactadental.com

*"WHEN I PLACE AN ORDER I CAN COUNT ON IT TO BE TO ME
WITHIN PROPER TIME AS PROMISED. THE CUSTOMER SER-
VICE IS OUTSTANDING. JUST REALLY LOVE THIS COMPANY!"*

- Sue Whitmer, RDA

our **focuses**

Provisionals

Restorative

Impressions

Composites

Provisional – beauty

When patients leave your office, you want them to look terrific whether they are wearing provisionals or final restorations. Our products will provide the aesthetics your patients want with the ease of use and cost savings that you deserve.

Restorative – strength

From core build up to fiber posts and our cements, great care has been taken to design products that not only benefit you, but your patients as well. We have combined technology and clinical experience to produce products that will surpass your expectations.

Impressions – accuracy

Find out how our nine viscosities of impression material can provide you with the flexibility for your technique while providing incredible accuracy for your lab. Impression material competitors can be gimmicky, but when you just want a material to work reliably, FRESH Bold Impression Material is your answer.

Composites – technology

There are many choices of composite in the dental market. Take a moment to review Nanoceram; a composite line that doesn't make you compromise. Using the latest technology, proven by independent research, we are able to offer you superior aesthetics, handling, longevity, and value.

EXACTA Temp **XTRA**

This is our best selling product! Our state of the art bis-acrylic material is available in 7 shades with a fast setting time of 1:20. Not only will you get life-like aesthetics but you can add EXACTA Temp to itself for repairs. Each kit contains 1 82g cartridge and 10 mixing tips.



EXACTA Cement **XTRA**

Zinc-oxide, non-eugenol formula will provide just the right amount of retention for your provisionals. Our formula will not discolor over time, gels nicely for cleanup, and removes cleanly from the prep. Available in an 11g handheld syringe or a 120g cartridge for dispensing guns.



EXACTA Flow Ic

For provisional repairs or general flowable composite needs, look no further than EXACTA Flow. It is shade matched perfectly to EXACTA Temp, so your repairs will be quick and easy. Cures in as little as 3 seconds.



Ordering Information

EXACTA Temp Xtra

SHADES: Bleach, B1, A1, A2, A3, A3.5, C2

CONTAINS: 1 x 82g cartridge and 10 mixing tips

PRICE: 1-3 kits - \$116.95, 4-9 kits - \$105.95, 10+ kits - \$100.95

EXACTA Cement Xtra

SYRINGE:

CONTAINS: 1 x 11g syringe and 10 mixing tips

PRICE: 1-3 kits - \$31.95, 4+ kits - \$28.95

CARTRIDGE:

CONTAINS: 1 x 120g cartridge and 25 mixing tips

PRICE: 1-3 kits - \$137.95, 4+ - \$127.95

EXACTA Flow Ic

SHADES: Bleach, B1, A1, A2, A3, A3.5, C2

CONTAINS: 1 x 2g syringe and 5 tips

PRICE: - \$12.95



FRESH BOLD VPS Impression Material

The single most important characteristic of an impression material is accuracy. The best test of that is how your crowns fit when seated; the rest is just marketing.

Big brands are now charging over - \$30 per cartridge with wild claims and citing 'study' after 'study'. The bottom line is if you want accurate impressions at a great price, FRESH BOLD VPS is for you.

Plan an upcoming case with FRESH Bold Impression Material and decide for yourself if paying twice as much for a brand name is worth it.

Final Impressions

FRESH BOLD VPS Impression Material provides you with the right materials for your method and provides accurate results consistently. With nine different viscosities to choose from, FRESH BOLD will help you capture the details you need to lead to a great restorations



Pre-Op Impressions

Why mess with alginate? Fresh BOLD Pre-Op Gel is a perfect material for a provisional matrix, bleaching trays, study models and other impression needs. This is a great alternative to alginate and is handy to have in emergency kits so you don't waste time mixing.



Bite Registrations

Whether you are looking for a traditional fluffy and quick bite registration, or rock-solid hardness, we have the choices for you. FRESH BOLD Bite Registration is our most popular and available in a 30 second or 2 minute version. Blu-Max has a rigid 93 Shore A hardness with a 60 second set time.



Ordering Information

FINAL IMPRESSIONS

Kits: (2 x 50ml cartridges and 6 mixing tips) - \$31.95

Bulks: (24 x 50ml cartridges and 72 mixing tips) - \$317.95

Heavy Body

Liquid Putty

Clear Bite Reg

Medium Body

Monophase

Light Body

Fast Set Only

Fast Set Only

Regular Set Only

Fast Set Only

Fast or Regular Set

Fast Set Only

Putty: (2 x 450ml jars) - \$105.95

Heavy Body: (380ml for mixing machine) - \$105.95

PRE-OP and BITE REGISTRATIONS

Kits: (2 x 50ml cartridges and 6 mixing tips) - \$25.95

Bulks: (24 x 50ml cartridges and 72 mixing tips) - \$254.95

Pre-Op Gel

Bite Registration

Fast or Regular Set

Fast or Regular Set

Blu-Max

Kits: (2 x 50ml cartridges and 6 mixing tips) - \$31.95

Bulks: (24 x 50ml cartridges and 72 mixing tips) - \$317.95

EXACTA Trays

"Rated excellent or good & purchase recommended by 100% of Evaluators"

These words changed our company. EXACTA Trays were our first product and once they were rated by the leading product evaluation organization, our credibility was established. The trays have simple yet functional design features: wide occlusal plane, rigid nylon-composite body, and each tray is individually wrapped.

Available in: Posterior (50), Anterior (36), Quadrant (40) , Full Arch (30), and Posterior Sideless (50).

Also available in an assorted box.

PRICE: 1-2 Boxes - \$42.95, 3-5 Boxes - \$39.95, 6+ Boxes - \$35.95



EXACTA Core-mft

A stable foundation is key to successful building. Your restorations will have just that. EXACTA Core-mft contains micro-fiber technology that provides great strength and physical characteristics of dentin. Dentists often remark on how cutting EXACTA Core-mft feels just like dentin. Flouride will help eliminate secondary carries and the material is very radiopaque to aid in diagnosis. Available in four shades: A2, Blue, White, and Gold and two dispensing systems.

Cartridge Kit: 1 x 25ml cartridge and 25 mixing tips & IO tips - \$95.95, 4+ kits - \$90.95

Syringe Kit: 4 x 4ml syringe and 20 mixing tips -- \$95.95, 4+ kits - \$90.95



Simplique

We decided early on that our permanent cement had to be simple and reliable. This is a self-adhesive resin cement that eliminates the need for etching and bonding for the cementation of almost all types of restorations. Simplique will reduce sensitivity and chairside time, and is a great choice for composite and all-ceramic restorations in the aesthetic zone. Its high bond strength and low expansion will provide great clinical results.

Starter Kit: 5 x 4ml syringes (3-A2, 1-White, 1-Translucent), Silane, and mixing tips - \$264.95

Refill kit: 1x 4ml syringe and mixing tips - \$63.95





Nanoceram Composite System

Sometimes you have to search high and low to find a gem. In this case, we found it in Greece. Nanoceram is manufactured by DMP Dental, a 27 year old company, and has been sold worldwide for years, but never had a distributor in the United States. Now, Nanoceram Universal, Flowable, and Posterior are available from EXACTA Dental Direct. This is truly a world-class material.

Often overlooked in the composite debate are two key factors - water sorbtion and solubility. The Dental Advisor did a research report on the subject and we were thrilled with the results. Nanoceram had lower sorbtion and solubility rates than many other well known brands. This means better and longer lasting restorations.

Beyond the science, this composite offer a beautiful result and maintains its glossy finish longer. Not too sticky, it is easily placed and sculpted and polishes quickly. Available in over 20 shades and in syringes or unit doses, Nanoceram is truly a practice changer.

Universal

Unit Dose: 20 x .25g - \$52.95

Syringe: 1 x 4g - \$48.95

Flowable

Syringe: 2 x 2ml - \$32.95

Posterior

Unit Dose: 20 x .25g - \$52.95

Syringe: 1 x 4ml - \$52.95



Bluespray

For CAD/CAM users who need to coat a surface for digital impressions, Bluespray provides visual definition with bright and crisp scans. Easy to apply and clean.

50ml can - \$79.95



StoneBite Scan

Rock hard bite registration pre-mixed with titanium dioxide for rigid, accurate, and easy to scan bite registrations.

Kit: 2 x 50ml cartridges, 12 mixing tips & ribbon tips - \$42.95



Opticide 3

No more measuring, mixing, or expiration. Opticide 3 has a three minute comprehensive kill time. Safe to use on all dental surfaces, Opticide will be a 'must' for your office. The wipes are large, blue surgical grade cloths that are perfect for your surface disinfection protocol.

100 Wipes - \$10.95

24oz. Spray Bottle - \$10.95

1 Gallon - \$26.95




Dispensing Guns

We offer three kinds of dispensing guns. The first is a 1:1 ratio, 50ml dispenser for impression materials and other 50ml cartridges. The second is a 4:1 ratio, 50ml dispenser for EXACTA Temp Xtra. Last, is a 1:1 ratio, 25ml dispenser for EXACTA Core cartridges. Our dispensing guns will work with other brands. Don't spend twice as much through a dealer!



Mixing Tips

We only distribute the authentic Sulzer-Mixpac dispensing system and carry a variety of different tips. Remember that mixing tips are universal and you don't need to buy the exact same tip for more from anywhere else.



Dispensing Gun Holder

If you want to cut down on chairside time, get all your guns in a row prior to the procedure. We also include a tip caddy that fits in the holder for even more organization.

Dispensing gun holder: \$63.95



Ordering Information

Dispensing Guns - \$84.95

25ml 1:1 Dispensing Gun

50ml 1:1 Dispensing Gun

50ml 4:1 Dispensing Gun

Mixing Tips

Impression Material:

48 tips - \$34.95, 6+ \$32.95

X-Large Purple

Large Green

Medium Pink

Small Yellow

Restorative Material:

25 tips - \$31.95, buy 3 bags, get 1 FREE

EXACTA Temp Xtra Tips

EXACTA Cement Xtra Tips

EXACTA Core mft Tips (cartridge only)

Intraoral Tips - 100 tips \$40.95

EXACTA Flow Tips - 50 tips \$10.95

Core Syringe Tips - 20 tips \$37.95

Mixing Machine Tips - 20 \$26.95

Notes:

EXACTA

Dental Direct

Buy Direct and Save | (800) 474-7665 | www.ExactaDental.com | info@exactadental.com