



Biolon





Biolon

Thermoforming blank made of PETG (Polyethylene terephthalate glycol) featuring an insulated film for high transparency. Tough elastic material that can be built up with acrylate. Suitable for aligners, bite and occlusion splints, as well as orthodontic appliances.

Application	Product	Thickness				
Aligner	Biolon	0.5 mm	0.625 mm	0.75 mm		
Orthodontic devices	Biolon	2 mm	3 mm			
Mouthguards	Biolon	2 mm				
Positioner / Retainer	Biolon	0.75 mm				
Production of temporaries	Biolon	0.5 mm	0.75 mm			
X-ray splints	Biolon	1 mm	1.5 mm	2 mm		
Splint Therapy	Biolon	0.75 mm	1 mm	1.5 mm	2 mm	3 mm