



# Drufosoft® color

Thermoforming blank made of ethyl vinyl acetate. Soft elastic with high recovery from deformation with a final hardness of 82 Shore A. Can be laminated and is ideally suited for individual sport mouthguard such as Drevo mouthguards, positioners, and duplication forms.

$\frac{1}{8}$	Thickness	Color	Pieces	Item no.
3 mm		Neon red	3	D42483X2
			10	D4248X2
3 mm		Neon green	3	D42483X3
			10	D4248X3
3 mm		Neon yellow	3	D42483X4
			10	D4248X4
3 mm		Blue	3	D42483X5
			10	D4248X5
3 mm		Yellow	3	D42483X6
			10	D4248X6
3 mm		Black	3	D42483X7
			10	D4248X7
3 mm		White	3	D42483X8
			10	D4248X8
3 mm		Pink	3	D42483X9
			10	D4248X9
3 mm		Red	3	D42483X17
			10	D4248X17
3 mm		Turquoise	3	D42483X18
			10	D4248X18
3 mm		Gold	3	D42483X19
			10	D4248X19
3 mm		Silver	3	D42483X20
			10	D4248X20
3 mm		Blue opaque	3	D42483X21
			10	D4248X21
3 mm		Green	3	D42483X22
			10	D4248X22
3 mm		Purple	3	D42483X23
			10	D4248X23
3 mm		Light blue	3	D42483X24
			10	D4248X24
3 mm		Bordeaux	3	D42483X25
			10	D4248X25
3 mm		Orange	3	D42483X26
			10	D4248X26