



## Drufosoft<sup>®</sup> transparent

Thermoforming blank made of ethyl vinyl acetate. Soft elastic with high recovery from deformation with a final hardness of 82 Shore A. Can be laminated and is ideally suited for individual sport mouthguard such as Dreve mouthguards, positioners, and duplication forms.

Thickness	Pieces	Item no.
1 mm	20	D4245
1 mm	50	D42451
1.5 mm	20	D4252
1.5 mm	50	D42531
2 mm	15	D4250
2 mm	50	D42511
3 mm	10	D4248
3 mm	30	D42481
4 mm	15	D42472
4 mm	75	D424901
5 mm	15	D4225