

# Thermofforming Blanks



Application	Product	Thickness				
Aligner	Biolon	0.5 mm	0.625 mm	0.75 mm		
Base plates	Drufoplast	2 mm				
	Drufoplast incl. tray handle	3 mm				
Bleaching	BioBleach® hard	0.5 mm				
	BioBleach® soft	1 mm				
	Drufosoft®	1 mm	1.5 mm			
Matrix	Drufolen H	0.6 mm				
	Drufolen W	0.45 mm	0.65 mm	0.8 mm		
Custom made trays	Drufoplast incl. tray handle	3 mm				
Orthodontic devices	Biolon	2 mm	3 mm			
Mouthguards	Biolon	2 mm				
	Drufosoft® transparent	2 mm	3 mm	4 mm		
	Drufosoft® color	3 mm				
Space maintainers	Drufolen W	0.65 mm	0.8 mm			
Positioner / Retainer	Biolon	0.75 mm				
	Drufosoft® transparent	1.5 mm	2 mm	3 mm	4 mm	5 mm
	Drufosoft® color	3 mm				
	Drufosoft® pro	3 mm				
Production of temporaries	Biolon	0.5 mm	0.75 mm			
	Drufolen H	0.6 mm	0.8 mm			
X-ray splints	Biolon	1 mm	1.5 mm	2 mm		
	Drufoplast inc. tray handle	3 mm				
Splint therapy	Biolon	0.75 mm	1 mm	1.5 mm	2 mm	3 mm
	Drufosoft® pro	3 mm				
	Kombiplast	2 mm	3 mm	4 mm		

